



Material	331562U
Material description	Buffer pH 9.22 (20°C) ± 0.02, tablets
CAS Number	1303-96-4
Additional information	

Characteristics	Specifications
Appearance:	.
White flat bevelled edge tablets	Passes test
Average tablet weight	979 - 1081 mg
Diameter	12.5 - 13.5 mm
Hardness	≥ 3 kp
pH (1 tablet; 100 ml)	9.20 - 9.24
Weight uniformity	927 to 1133 mg

Signature	Additional information
We certify that this batch conforms to the specifications listed above. This document has been produced electronically and is valid without a signature. VWR International	Weight uniformity: max. 10 % outside 953 - 1107 mg; none outside 927 - 1133 mg