



Specification

1.01187.1000 Ammonium nitrate EMPLURA®

Spec. Values		
Assay (alkalimetric)	≥ 95.0	%
Identity	passes test	
Chloride (Cl)	≤ 0.002	%
Phosphate (PO ₄)	≤ 0.0005	%
Sulfate (SO ₄)	≤ 0.01	%
Heavy metals (as Pb)	≤ 0.001	%
Ca (Calcium)	≤ 0.005	%
Fe (Iron)	≤ 0.0005	%
Mg (Magnesium)	≤ 0.001	%
Residue on ignition (as sulfate)	≤ 0.05	%

Claudia Wiegand

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.