



CERTIFICATE OF ANALYSIS

DESCALTM

Lot # 9J110

Manufactured October 2019

Recommended Retest date is two years from the date of manufacture

SIZE	1L (34oz)	3.79 L(1gal)	20L (5.28gal)	115L (30gal)
DECON #	4001	4004	4002	4030

TEST	SPECIFICATION	TEST RESULT
Appearance	Clear and colorless	Clear Liquid
Specific Gravity	1.165 - 1.185	1.175
pH	1.0 - 2.0	1.5

Reviewed and Approved by:

Date: October 2019

A handwritten signature in blue ink that reads 'Karen Wend'.

Karen Wend
Manager, Quality and Technical Service

BSE/TSE Statement: There are no materials of animal origin in the raw materials or in the process aides used to manufacture any Decon products.