



KIMTECH SCIENCE®
KIMWIPES® Ex-L Delicate Task Wipers
Certificate of Conformance

This document certifies that KIMTECH SCIENCE® KIMWIPES® Ex-L Delicate Task Wipers conform to Kimberly-Clark Corporation's internal specifications for this product as stated below. We use a system of in-process and audit inspections to assure that only product which meets specification is released for shipment to customers.

Certified by: *Mpekele White*

Comments: MO4360

Code	Type	Color/ Colour	Material	Size	Count
34120	Pop-Up Box	White	1 ply tissue	8.2" x 4.4" 20.8 x 11.2 cm	286 wipers/box 30 boxes/case
34133	Pop-Up Box	White	1 ply tissue	12.3" x 11.22" 31.2 x 28.5 cm	198 wipers/box 15 boxes/ case
34155	Pop-Up Box	White	1 ply tissue	8.2" x 4.4" 20.8 x 11.2 cm	286 wipers/box 60 boxes/case
34256	Pop-Up Box	White	1 ply tissue	16.4" x 14.43" 41.7 x 36.7 cm	144 wipers/box 15 boxes/case

Physical property data is as produced, subsequent storage conditions and the age of the product, which are beyond Kimberly-Clark's control, may affect the properties of the product as received.

	Average	Test Method
CD Dry Tensile Strength		
g/3" strip	1,500	STP-578
CD Wet Tensile Strength		
g/3" strip	425	STP-578
MD Dry Tensile Strength		
g/3" strip	1500	STP-578
Water Capacity (1 ply)		
g/4" x 4"	1.1	RTM-6430
Oil Capacity		
g/4" x 4"	0.60	RTM-6430
Gelbo Lint		
#of particulates >50 micrometers	35	RTM-6176
Silicone		
ppm	15	BF3-GC

Should you require additional information, please contact our Technical Help Line at 888/346-4652.