



Product Specification Sheet

VWR SQUARE GLASS CUVETTES, LIGHT PATH 100 mm

VWR Cat.No.	VWRI634-6017
Product Name	SQUARE.GLASS CUVETTE (100mm)
Serial Number	634-6017
Material	B270 Optical Glass
Material Thickness	1.25±0.1 (mm)
Size	102.5*12.5*45 L*W*H (mm)
Optical Path Length	100±0.15 (mm)
Optical Path Width	10±0.1 (mm)
Cubic Measure	35 (ml)
Transmittance	≥79% (at 700nm)
Package	1pcs/box

VWR International bvba/sprl, Haasrode Research Park Zone 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>

Date of issue: Dec.03rd 2012