

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR® Ice Pans with or without lid, Stackable, Leakproof
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2022-03-04

Without Lid

Catalogue Number	Capacity [L]	Colour	Dimensions (W×D×H) [mm]
216-1525	1	Black	215 × 188 × 118
216-1526	1	Blue	215 × 188 × 118
216-1527	1	Green	215 × 188 × 118
216-1528	1	Red	215 × 188 × 118
216-1529	4	Black	265 × 380 × 156
216-1530	4	Blue	265 × 380 × 156
216-1531	4	Green	265 × 380 × 156
216-1532	4	Red	265 × 380 × 156
216-1533	9	Black	356 × 503 × 155
216-1537	9	Blue	356 × 503 × 155
216-1538	9	Green	356 × 503 × 155
216-1539	9	Red	356 × 503 × 155

With Lid

Catalogue Number	Capacity [L]	Colour	Dimensions (W×D×H) [mm]
216-1493	1	Black	215 × 188 × 118
216-1494	1	Blue	215 × 188 × 118
216-1495	1	Green	215 × 188 × 118
216-1496	1	Red	215 × 188 × 118
216-1497	4	Black	265 × 380 × 156

Catalogue Number	Capacity [l]	Colour	Dimensions (W×D×H) [mm]
216-1498	4	Blue	265 × 380 × 156
216-1499	4	Green	265 × 380 × 156
216-1503	4	Red	265 × 380 × 156
216-1506	9	Black	356 × 503 × 155
216-1507	9	Blue	356 × 503 × 155
216-1508	9	Green	356 × 503 × 155
216-1509	9	Red	356 × 503 × 155

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Test Name	Test Reference	Accuracy	Result
Leak Test	QWI-QC-06	General Insp. Level- I, AQL- 1.5	Leak proof

Product Specifications

The product is in conformance with the requirements of the following standards and regulations:

Material	Polyurethane.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).

Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The above products are not made from natural rubber latex.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Room temperature or ambient storage condition.
Shelf Life	Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature