



VWR® BEAKERS WITH HANDLE, POLYPROPYLENE

As per the manufacturer, the below product meets the following criteria:

VWR Catalog Number	75800-262
Lot Number	V17-011119
Description	VWR BEAKER AUTOCLAVABLE W/HDL PP 5000ML
Date of Manufacture	November 01, 2019
Country of Origin	India.

Quality System Compliance

Products are manufactured, inspected and found to be in compliance with product and quality specification requirements as documented in the ISO 9001:2015 & ISO13485:2016 Quality Management System. Inspection records are reviewed and signed for product release.

QC Testing

Representative production samples are collected and inspected in accordance with current applicable standard protocols.

Test Name	Test Reference	Acceptable Quality Level	Result
Autoclaving	QWI-QC-07	Special Insp. Level-S-1, AQL- 1.5	Autoclavable
Dimension	QWI-QC-23	Special Insp. Level-S-1, AQL- 1.5	Within tolerance

Product Specifications

Capacity (mL): 5000 **Graduations (mL):** 1000 **Sub divisions:** 500 m **Weight:** 417 ± 20 g

Tolerance: ± 250 mL **Maximum Top Outer Diameter:** 260 ± 2 mm **Height:** 262 ± 2 mm

Material: Polypropylene (PP), 100 % Virgin, FDA 21 CFR Compliance, Food grade.

BSE/TSE Statement: The above-mentioned product is manufactured without the use of Raw materials of animal origin. Therefore, BSE/TSE issue does not concern this product to the best of the manufacturer's knowledge.

Non-Pyrogenic: Not Applicable.

DNase & RNase Free: Not Applicable.

Cytotoxicity: Not Applicable.

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Storage Conditions: Can be stored in room temperature in normal warehouse conditions.

Signed:



Jamie Ethier
VP Global Quality
VWR, Part of Avantor

Date: October 13, 2020