

Certificate of Conformity

VWR Catalogue Number	113-0716
Description	Face masks, 3 ply, non woven PP, with earloops, blue
Lot Number	D9210097
Date of Manufacture (yyyy-mm-dd)	2021-12-25
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2024-02-23

Quality System Compliance

Products are manufactured under the **ISO 9001:2008 & ISO 13485:2003** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

Bacteria filtration rate (%)

SPEC: >90%

Test Result: >90%

Material: SPP 25g/m2 White Inner PLY

	Test Method	MD	CD
Break Tensile Properties [N]	ASTM D5034	69.40	64.77
Elongation rate [%]	ASTM D5034	122.32	117.37

Material: SPP 20g/m2 Light Blue Outer PLY

	Test Method	MD	CD
Break Tensile Properties [N]	ASTM D5034	30.28	25.47
Elongation rate [%]	ASTM D5034	106.34	104.55

Product Specifications

Material	Non-woven fabric 20g PP/20g MB/25g PP.
Density	Front layer: 20 g/m ² Back layer: 25 g/m ² Centre: 20 g/m ²
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Room temperature, dry and shady.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.