




TYPE I FLINT

NOTE

VIAL MAY BE USED
FOR FREEZE DRYING

① ADDED 8-16-88

CUSTOMER:	ALCAN PACKAGING  Tubing Americas 1200 North Tenth Street Millville, N.J. 08332-2032 USA				DRAWN BY: EM	DATE DRAWN 11-7-85
TITLE: SERUM VIAL					SCALE: FULL OR AS NOTED	71000850
FINISH: 20mm A/S w/BLOWBACK					REDRAWN 3-24-94	DEV.#
CAPACITY SPECIFICATION:	APPROVAL	ENGINEERING	MANUFACTURING	QUALITY ASSURANCE		
PRACTICAL FILL: 5CC	INITIAL & DATE	CEP	10-23-03	DGB	10-26-03	LR 10-23-03
				ARTICLE # 2905-B8B		