







Tolerances unless otherwise specified:
Without decimal: ±0.5 mm
1 place decimal: ±0.1 mm
2 place decimal: ±0.05 mm

Undefined radius: R 0.20 mm

Roskilde Site · Kamstrupvej 90 · P.O.Box 280
DK-4000 Roskilde · Denmark
Tel +45 4631 2000 · Fax +45 4631 2042
E-mail: info.nunc@thermofisher.com

www.nuncbrand.com

SUBJECT:		DRAWN BY/DATE:		CHANGE ORDER NO.:		VALID FROM:		APPROVED BY/DATE:		FORMAT:
Cryo tube 1.0 ml. w. internal thread and ribs		PS/23.03.93		CO-1181		2009.03.26		JR 2009.03.12		A3
ITEM NUMBER:		MATERIAL:		DRAWING NO.:		VERSION:		PAGE TYPE:		SCALE:
-		PP		2028		10-a		-C		5:1, 1:1
										METRIC
										First Angle