




Dimensions (Unit: mm)



VWR® 5mL Micro Tube

Description: 5mL Micro Tube, RNase/DNase & Pyrogen Free

Purpose: Used for sample storage, operation and centrifugation

Materials

Tube/Cap: PP (Polypropylene)

Features

- Made from high grade gamma resistance polypropylene
- Autoclavable at 121°C and freezable to -80°C
- Engraved graduation ensure accuracy
- Flat and frost caps surface together with smooth and frosted body surface
- Provide easy and legible mark area
- Tubes can be centrifuged to 25,000g RCF
- DNase/RNase-free and Non-pyrogenic
- Sterilized by irradiation SAL10⁻⁶ (ISO11137)
- Shelf life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® 5mL Micro Tube

VWR EU Cat. No.	Capacity	Tube Color	Bottom Style	Graduation	Sterilization	Packaging
525-0666	5.0 mL	Natural	Conical	Yes	Yes	20/PK, 200/case
525-0660	5.0 mL	Natural	Conical	Yes	No	200/PK
525-0661	5.0 mL	Blue	Conical	Yes	No	200/PK
525-0662	5.0 mL	Yellow	Conical	Yes	No	200/PK
525-0663	5.0 mL	Green	Conical	Yes	No	200/PK
525-0664	5.0 mL	Rose	Conical	Yes	No	200/PK
525-0665	5.0 mL	Black	Conical	Yes	No	200/PK