

Material	24611.366
Material description	n-Hexane
Grade	TECHNICAL
CAS Number	110-54-3
Molecular formula	$\text{H}_3\text{C}(\text{CH}_2)_4\text{CH}_3$
Molecular mass	86.18

Characteristics	Specifications
Assay	$\geq 95 \%$
Boiling point	68 - 70 °C
Benzene	$\leq 100 \text{ ppm}$
Evaporation residue	$\leq 20 \text{ ppm}$

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.