

Certificate of Quality

VWR Catalogue Number	82027-584
Description	VWR® Dissecting Scissors, Curved Tip, 4 ¹ / ₂ " (114 mm)
Country of Origin	Manufactured in Pakistan
Date of Issue (yyyy-mm-dd)	2022-09-02

Quality System Compliance

Products are manufactured under the [ISO 9001:2008 & ISO13485:2003 & F-899](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Following tests performs as per procedure:

- Material Chemical Analysis 7-element Test.
- Boil Test
- Passivation

Product Specifications

Material	AISI 410 type.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.

Latex Statement	The product is Latex free.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Store in a dry place.
Shelf Life	No Shelf Life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.