



Product Specification Sheet

VWR CLEANROOM GOGGLES, STERILE, SINGLE USE, ANTI-FOG, ANTI-SCRATCH

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No: 111-0615

CERTIFICATION: EN 166:2001 Optical Class 1.

QUALITY SYSTEMS: Manufactured in a facility holding ISO 9001:2008 certification and packed in a NEBB certified ISO Class 4 (Class 10) cleanroom.

PACKAGING: One piece per Tyvek/PE bag; 10 pieces per outer Tyvek/PE bag; Six outer bags per carton liner; One carton liner per carton (60 Pieces).

COMPATIBILITY: ISO Class 4 and EU GMP Grade A/B critical environments.

CLASSIFICATION: Category 2 - Personal Protective Equipment conforming to PPE Directive 89/686/EEC.

OTHER CHARACTERISTICS:

- Sterile cleanroom goggles with a lightweight frame, indirect ventilation and a non linting latex-free head band.
- Toughened polycarbonate lens with anti-fog and anti-scratch coating.
- Lightweight PVC body with silicone rubber strap provide comfort.
- Ethylene Oxide (EtO) sterilised, fully validated to ISO 11135-1-2007.

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **June 2013**