



CERTIFICATE OF MANUFACTURE/ANALYSIS

Product Size:	1 Gallon	Product Name:	LpH III se Phenolic Disinfectant
Manufacture Date:	April 28, 2022	Product ID:	6470-08
Exp. Date:	April 1, 2024	Lot No.:	318767

I. The product listed above is made per the current issue of the following documentation:

Raw Material Quality Control Specifications
Finished Product Master Formula
Manufacturing Instructions
Packaging Specifications
Finished Product Quality Control Specifications

Copies of this documentation are on file at STERIS Corporation, St. Louis Operations and will be provided as required to the applicable regulatory agencies.

II. The product listed above has been analyzed. The release specifications and results are shown below.

<u>Test Name</u>	<u>Method No.</u>	<u>Results</u>	<u>Specification</u>
Color	110	Light Straw	Light Straw to Dark Amber
Clarity	107	Clear	Clear to Slightly Hazy
Clarity (Nephelometer)	225	5	≤ 20
pH, 1% w/w	106	2.1	1.6 - 2.5
Cloud Point	128	< 15°C	<15°C
OPP Assay	LC52	6.10%	5.82 - 6.30%
OBPCP Assay	LC52	13.44%	12.67 - 13.65%
Specific Gravity (@ 25°C)	114	1.116	1.108 – 1.127

Based on these results the product is released for shipment to all STERIS Corporation distribution points.

Signed: *Erison*

Date: May 2, 2022

Quality Representative