

DWK Life Sciences LLC.
1501 North Tenth Street
Millville NJ 08332

Certificate# 2100645739-01

Sold to:
VWR INTERNATIONAL LLC
P.O. BOX 2617
SECAUCUS NJ 07094

Ship to:
VWR VISALIA RDC
8711 W. RIGGINS AVE
VISALIA CA 93291

Certificate of Compliance:

In accordance with applicable ISO 9001 standards, cGMP's, and Company requirements, DWK Life Sciences LLC. hereby certifies the manufactured item described above (collectively, the "Product") that:

1. General: Product has been given appropriate visual, and/or mechanical and/or electronic inspection and the results of these inspections indicate that the product meets all agreed upon contractual specifications at the time of its shipment.
2. Glass: We certify that the product made from:
TYPE I Borosilicate Glass meets the requirements for the glass grain test limits for Type I Glass when tested according to the current revisions of the USP < 660 > and EP3.2.1 Pharmacopeias
Arsenic Extraction this ware will meet the requirements as described in the current revisions of USP < 660 > and EP3.2.1 Pharmacopeias for Arsenic Extraction.
TYPE III Soda Lime Glass meets the requirements for the glass grain test limits according to the current revisions of the USP < 660 > and EP3.2.1 Pharmacopeias.
Amber Glass conforms to the light transmission limits specified in the current revision of the USP and EP Pharmacopeias.
Hydrolytic Resistance when indicated above this ware conforms to the current revisions USP < 660 > Surface Test and EP 3.2.1.A - Test for surface hydrolytic resistance.
Heavy Metals are not intentionally added to the glass compositions. Thus Contents of lead, cadmium, mercury and hexavalent chrome are in compliance with the limit values of 100 parts per million by weight for US (TPCH and CONEG) and EC regulations.
Transmissible and Bovine Spongiform Encephalopathy (TSE/BSE) - We certify that containers made from glass are free of any nonconforming animal derivatives.
3. Plastic: We certify that product made from resin material has been manufactured under controlled conditions and is in compliance with all applicable procedures and specifications.
4. Rubber: We certify that product made from rubber has been manufactured under controlled conditions and is in compliance with all applicable procedures and specifications.
5. Aluminum: We certify that product made from aluminum has been manufactured under controlled conditions and is in compliance with all applicable procedures and specifications.

Certified by: Elizabeth Wonderlin

Date: 2023-09-05 10:33

Title: QC Technician



DWK Life Sciences LLC.
1501 North Tenth Street
Millville NJ 08332

Certificate of Compliance

Certificate# 2100645739-01

Customer Material Number	: 89000-940	Customer Purchase Order	: 4514415537
DWK Part Number	: 9143950008		
DWK Description	: CAP, GL45, BLU PP		
Cert Date	: 2023-09-05	Shipment Quantity	: 5 CAS
EPHR	: No		

Batch/Lot Number	: 08/21/2023	Batch/Lot Quantity	: 5 CAS
------------------	--------------	--------------------	---------

Customer Material Number	: 89000-940	Customer Purchase Order	: 4514426038
DWK Part Number	: 9143950008		
DWK Description	: CAP, GL45, BLU PP		
Cert Date	: 2023-09-05	Shipment Quantity	: 5 CAS
EPHR	: No		

Batch/Lot Number	: 08/21/2023	Batch/Lot Quantity	: 5 CAS
------------------	--------------	--------------------	---------

Customer Material Number	: 89000-940	Customer Purchase Order	: 4514442443
DWK Part Number	: 9143950008		
DWK Description	: CAP, GL45, BLU PP		
Cert Date	: 2023-09-05	Shipment Quantity	: 15 CAS
EPHR	: No		

Batch/Lot Number	: 08/21/2023	Batch/Lot Quantity	: 15 CAS
------------------	--------------	--------------------	----------

Customer Material Number	: 89000-940	Customer Purchase Order	: 4514565408
DWK Part Number	: 9143950008		
DWK Description	: CAP, GL45, BLU PP		
Cert Date	: 2023-09-05	Shipment Quantity	: 11 CAS
EPHR	: No		

Batch/Lot Number	: 08/21/2023	Batch/Lot Quantity	: 11 CAS
------------------	--------------	--------------------	----------

Customer Material Number	: 89000-940	Customer Purchase Order	: 4514675463
DWK Part Number	: 9143950008		
DWK Description	: CAP, GL45, BLU PP		
Cert Date	: 2023-09-05	Shipment Quantity	: 15 CAS
EPHR	: No		



Certificate of Compliance

DWK Life Sciences LLC.
1501 North Tenth Street
Millville NJ 08332

Certificate# 2100645739-01

Batch/Lot Number	: 08/21/2023	Batch/Lot Quantity	: 15 CAS
------------------	--------------	--------------------	----------

Customer Material Number	: 89000-940	Customer Purchase Order	: 4514680976
DWK Part Number	: 9143950008		
DWK Description	: CAP, GL45, BLU PP		
Cert Date	: 2023-09-05	Shipment Quantity	: 22 CAS
EPHR	: No		

Batch/Lot Number	: 08/21/2023	Batch/Lot Quantity	: 22 CAS
------------------	--------------	--------------------	----------

Customer Material Number	: 89002-050	Customer Purchase Order	: 4515299868
DWK Part Number	: 9280170011		
DWK Description	: FLASK, VOLUME, CLASSA, MICRO, CYL, 5ML		
Cert Date	: 2023-09-05	Shipment Quantity	: 1 CAS
EPHR	: No		

Batch/Lot Number	: 220188184	Batch/Lot Quantity	: 1 CAS
------------------	-------------	--------------------	---------

Customer Material Number	: 89000-940	Customer Purchase Order	: 4515548246
DWK Part Number	: 9143950008		
DWK Description	: CAP, GL45, BLU PP		
Cert Date	: 2023-09-05	Shipment Quantity	: 8 CAS
EPHR	: No		

Batch/Lot Number	: 08/21/2023	Batch/Lot Quantity	: 8 CAS
------------------	--------------	--------------------	---------

Customer Material Number	: 89003-530	Customer Purchase Order	: 4515658062
DWK Part Number	: 9370340026		
DWK Description	: PIPET, SERO, COTTON, LARGE, WHT, 25X0.1ML		
Cert Date	: 2023-09-05	Shipment Quantity	: 2 CAS
EPHR	: No		

Batch/Lot Number	: 220187848	Batch/Lot Quantity	: 2 CAS
------------------	-------------	--------------------	---------

Customer Material Number	: 66015-020	Customer Purchase Order	: 4515774793
DWK Part Number	: 225534		
DWK Description	: 8 mL Btl, Sample, Type III Clr, Snap Ca		



Certificate of Compliance

DWK Life Sciences LLC.
1501 North Tenth Street
Millville NJ 08332

Certificate# 2100645739-01

Cert Date	: 2023-09-05	Shipment Quantity	: 5 CAS
EPHR	: No		

Batch/Lot Number	: 0000050103	Batch/Lot Quantity	: 5 CAS
------------------	--------------	--------------------	---------

Customer Material Number	: 16620-004	Customer Purchase Order	: 4515817005
DWK Part Number	: 9151100000		
DWK Description	: BOTTLE, SPEC GRAVITY, HUBBARD, 24ML		
Cert Date	: 2023-09-05	Shipment Quantity	: 1 CAS
EPHR	: No		

Batch/Lot Number	: 220188464	Batch/Lot Quantity	: 1 CAS
------------------	-------------	--------------------	---------

Customer Material Number	: 93000-706	Customer Purchase Order	: 4515822484
DWK Part Number	: 9721205002		
DWK Description	: PIPET, SERO, GLASS, BULK, 10X0.1ML, VWR		
Cert Date	: 2023-09-05	Shipment Quantity	: 3 CAS
EPHR	: No		

Batch/Lot Number	: 23T17M	Batch/Lot Quantity	: 3 CAS
------------------	----------	--------------------	---------

Customer Material Number	: 93000-706	Customer Purchase Order	: 4515891986
DWK Part Number	: 9721205002		
DWK Description	: PIPET, SERO, GLASS, BULK, 10X0.1ML, VWR		
Cert Date	: 2023-09-05	Shipment Quantity	: 7 CAS
EPHR	: No		

Batch/Lot Number	: 23T17M	Batch/Lot Quantity	: 7 CAS
------------------	----------	--------------------	---------

Customer Material Number	: 89524-188A	Customer Purchase Order	: 4515939037
DWK Part Number	: 89524-188A		
DWK Description	: VWR 2ML 9MM S/T PP VIAL 100/PK		
Cert Date	: 2023-09-05	Shipment Quantity	: 27 PAK
EPHR	: No		

Batch/Lot Number	: 4600088397	Batch/Lot Quantity	: 27 PAK
------------------	--------------	--------------------	----------

Customer Material Number	: 16150-394	Customer Purchase Order	: 4515938983
--------------------------	-------------	-------------------------	--------------



Certificate of Compliance

DWK Life Sciences LLC.
1501 North Tenth Street
Millville NJ 08332

Certificate# 2100645739-01

DWK Part Number	: 220946		
DWK Description	: 32 oz Btl, Type III Amb, Safety Ctd, Pe Lnr		
Cert Date	: 2023-09-05	Shipment Quantity	: 9 CAS
EPHR	: No		

Batch/Lot Number	: 0000050130	Batch/Lot Quantity	: 9 CAS
------------------	--------------	--------------------	---------

Customer Material Number	: 46610-720A	Customer Purchase Order	: 4515997606
DWK Part Number	: 46610-720A		
DWK Description	: VWR CLOSURE 11MM SN PTFE PK100		
Cert Date	: 2023-09-05	Shipment Quantity	: 9 CAS
EPHR	: No		

Batch/Lot Number	: 293468	Batch/Lot Quantity	: 9 CAS
------------------	----------	--------------------	---------

Customer Material Number	: 46610-700A	Customer Purchase Order	: 4515997606
DWK Part Number	: 46610-700A		
DWK Description	: VWR VIAL BLK PP PTFE/SIL CS100		
Cert Date	: 2023-09-05	Shipment Quantity	: 27 CAS
EPHR	: No		

Batch/Lot Number	: 298036	Batch/Lot Quantity	: 27 CAS
------------------	----------	--------------------	----------

Customer Material Number	: 46610-716A	Customer Purchase Order	: 4515997606
DWK Part Number	: 46610-716A		
DWK Description	: VWR SEAL 11MM PTFE/SILCN PK100		
Cert Date	: 2023-09-05	Shipment Quantity	: 2 CAS
EPHR	: No		

Batch/Lot Number	: 300105	Batch/Lot Quantity	: 2 CAS
------------------	----------	--------------------	---------

Customer Material Number	: 82028-434A	Customer Purchase Order	: 4515997606
DWK Part Number	: 82028-434A		
DWK Description	: VWR CAP 10-425 PP BLACK PK100		
Cert Date	: 2023-09-05	Shipment Quantity	: 11 PK
EPHR	: No		

Batch/Lot Number	: 300110	Batch/Lot Quantity	: 11 PK
------------------	----------	--------------------	---------



DWK Life Sciences LLC.
1501 North Tenth Street
Millville NJ 08332

Certificate of Compliance

Certificate# 2100645739-01

Customer Material Number	: 89239-018A	Customer Purchase Order	: 4515997606
DWK Part Number	: 89239-018A		
DWK Description	: VWR CAP BL PP PTFE/SIL/PTFE		
Cert Date	: 2023-09-05	Shipment Quantity	: 4 CAS
EPHR	: No		

Batch/Lot Number	: 297528	Batch/Lot Quantity	: 4 CAS
------------------	----------	--------------------	---------

Customer Material Number	: 89524-190A	Customer Purchase Order	: 4515997606
DWK Part Number	: 89524-190A		
DWK Description	: VWR 300UL 9MM S/T PP VIAL 100/PK		
Cert Date	: 2023-09-05	Shipment Quantity	: 11 PAK
EPHR	: No		

Batch/Lot Number	: 4600109779	Batch/Lot Quantity	: 11 PAK
------------------	--------------	--------------------	----------

Customer Material Number	: 76047-554A	Customer Purchase Order	: 4515997606
DWK Part Number	: 76047-554A		
DWK Description	: VWR CAP NL TS SEPTA RED 9MM PK 100		
Cert Date	: 2023-09-05	Shipment Quantity	: 20 PK
EPHR	: No		

Batch/Lot Number	: 298759	Batch/Lot Quantity	: 20 PK
------------------	----------	--------------------	---------

Customer Material Number	: 60824-433	Customer Purchase Order	: 4515997554
DWK Part Number	: 9450480011		
DWK Description	: TUBE, CULTURE, PR, 25X250MM		
Cert Date	: 2023-09-05	Shipment Quantity	: 1 CAS
EPHR	: No		

Batch/Lot Number	: 146347	Batch/Lot Quantity	: 1 CAS
------------------	----------	--------------------	---------

Customer Material Number	: 89000-808	Customer Purchase Order	: 4515997554
DWK Part Number	: 9282600004		
DWK Description	: CRUCIBLE, GOOCH, LOW FORM, FRITTED, M, 30ML		
Cert Date	: 2023-09-05	Shipment Quantity	: 1 CAS
EPHR	: No		



DWK Life Sciences LLC.
1501 North Tenth Street
Millville NJ 08332

Certificate of Compliance

Certificate# 2100645739-01

Batch/Lot Number	: 220184904	Batch/Lot Quantity	: 1 CAS
------------------	-------------	--------------------	---------

Customer Material Number	: 89174-488	Customer Purchase Order	: 4515997554
DWK Part Number	: W240860		
DWK Description	: CAPCLUSTER MATS FOR MICROTUBES,TPE,NAT		
Cert Date	: 2023-09-05	Shipment Quantity	: 2 CAS
EPHR	: No		

Batch/Lot Number	: 0000048314	Batch/Lot Quantity	: 2 CAS
------------------	--------------	--------------------	---------

Customer Material Number	: 66010-695	Customer Purchase Order	: 4515997554
DWK Part Number	: 9738089907		
DWK Description	: CLOSURE,PPL,PTFE/SILICONE		
Cert Date	: 2023-09-05	Shipment Quantity	: 1 CAS
EPHR	: No		

Batch/Lot Number	: 220184688	Batch/Lot Quantity	: 1 CAS
------------------	-------------	--------------------	---------

Customer Material Number	: 66022-525	Customer Purchase Order	: 4515997554
DWK Part Number	: 868806		
DWK Description	: 50 Pos Rack, PP, 28.1 mm Open Id		
Cert Date	: 2023-09-05	Shipment Quantity	: 4 CAS
EPHR	: No		

Batch/Lot Number	: 4600077816	Batch/Lot Quantity	: 4 CAS
------------------	--------------	--------------------	---------

Customer Material Number	: 97047-294	Customer Purchase Order	: 4515997554
DWK Part Number	: 11-1400		
DWK Description	: Crimp uLvial 12x32mm Clear		
Cert Date	: 2023-09-05	Shipment Quantity	: 4 PK
EPHR	: No		

Batch/Lot Number	: 0000049523	Batch/Lot Quantity	: 4 PK
------------------	--------------	--------------------	--------

Customer Material Number	: 97047-994	Customer Purchase Order	: 4515997554
DWK Part Number	: 20-2400		
DWK Description	: Vial Flat Btm 22x75mm Clr 20mL		



Certificate of Compliance

DWK Life Sciences LLC.
1501 North Tenth Street
Millville NJ 08332

Certificate# 2100645739-01

Cert Date : 2023-09-05
EPhR : No

Shipment Quantity : 12 PK

Batch/Lot Number : 0000049027

Batch/Lot Quantity : 6 PK

Batch/Lot Number : 0000049624

Batch/Lot Quantity : 6 PK

Customer Material Number : 16171-818

Customer Purchase Order : 4515997554

DWK Part Number : 224182-01

DWK Description : 13 mm Seal, Ctr Tear-Out, Alum Nat, UnIn

Cert Date : 2023-09-05

Shipment Quantity : 2 CAS

EPhR : No

Batch/Lot Number : 0000048473

Batch/Lot Quantity : 2 CAS

Customer Material Number : 25485-029

Customer Purchase Order : 4515997554

DWK Part Number : 374320

DWK Description : Sinkers For 3 mm Tubing Assy

Cert Date : 2023-09-05

Shipment Quantity : 15 CAS

EPhR : No

Batch/Lot Number : 0000050302

Batch/Lot Quantity : 15 CAS

Customer Material Number : 66009-455

Customer Purchase Order : 4515997554

DWK Part Number : 9609580003

DWK Description : VIAL, KG33, 30X125MM, 24-400, NO CAP, 60ML

Cert Date : 2023-09-05

Shipment Quantity : 4 CAS

EPhR : No

Batch/Lot Number : 2-1-9

Batch/Lot Quantity : 4 CAS

Customer Material Number : 97048-652

Customer Purchase Order : 4515997554

DWK Part Number : F5-AL10

DWK Description : PP Syringe 10mL Luer Lock

Cert Date : 2023-09-05

Shipment Quantity : 2 PK

EPhR : No



Certificate of Compliance

DWK Life Sciences LLC.
1501 North Tenth Street
Millville NJ 08332

Certificate# 2100645739-01

Batch/Lot Number	: 0000037243	Batch/Lot Quantity	: 2 PK
------------------	--------------	--------------------	--------

Customer Material Number	: 490009-886	Customer Purchase Order	: 4515997554
DWK Part Number	: 292392809		
DWK Description	: SCREWCAPS GL45 BLUE PP		
Cert Date	: 2023-09-05	Shipment Quantity	: 3 PK
EPHR	: No		

Customer Material Number	: 89003-570	Customer Purchase Order	: 4515997554
DWK Part Number	: 9737500001		
DWK Description	: TUBE,DCT,NO CAP,15-415,16X100MM		
Cert Date	: 2023-09-05	Shipment Quantity	: 4 CAS
EPHR	: No		

Batch/Lot Number	: 220184735	Batch/Lot Quantity	: 4 CAS
------------------	-------------	--------------------	---------

Customer Material Number	: 93000-680	Customer Purchase Order	: 4515997554
DWK Part Number	: 9721005002		
DWK Description	: PIPET,SERO,GLASS,MULTI,10X0.1ML,VWR		
Cert Date	: 2023-09-05	Shipment Quantity	: 6 CAS
EPHR	: No		

Batch/Lot Number	: CC69G	Batch/Lot Quantity	: 6 CAS
------------------	---------	--------------------	---------

Customer Material Number	: KT294000-0010	Customer Purchase Order	: 4515997554
DWK Part Number	: 9294000002A		
DWK Description	: FLK RB 1N 10ML 14/20 BW		
Cert Date	: 2023-09-05	Shipment Quantity	: 3 CAS
EPHR	: No		

Batch/Lot Number	: 31/10/20	Batch/Lot Quantity	: 1 CAS
Batch/Lot Number	: 4600051498	Batch/Lot Quantity	: 2 CAS