

Product Specification

Test	Specification
The following results are based on the bulk bonded phase:	
Carbon (C)(%)	For Information Only
Nitrogen (N)(%)	For Information Only
Surface Coverage (ug/m2), Based on %C	350 – 450
Surface Coverage (ug/m2), Based on %N	85 – 105
The following results are based on the silica gel raw material:	
Average Particle Diameter, um (APD) – via Malvern (100 mm lens)	47 – 61
The following results are based on the solid phase extraction columns:	
Extraction Residue	<= 0.2 %

For Laboratory, Research or Manufacturing use.
6mL Solid Phase Extraction Columns
500 mg per Column



Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000