

Material	26846.292
Material description	Potassium iodide
Grade	AnalaR NORMAPUR
CAS Number	7681-11-0
Molecular formula	KI
Molecular mass	166.00

Characteristics	Specifications
Assay	≥ 99.5 %
Solution in water	Passes test
Alkalinity	≤ 0.002 meq/g
pH (20°C; 5 %)	5.7 - 9.2
Heavy metals (as Pb)	≤ 2 ppm
Insolubility in water	≤ 50 ppm
Loss on drying (150°C)	≤ 0.2 %
Reducing substances (as I)	≤ 10 ppm
Total N (Nitrogen)	≤ 10 ppm
Cl + Br (as Cl)	≤ 100 ppm
IO ₃ (Iodate)	≤ 2 ppm
PO ₄ (Phosphate)	≤ 10 ppm
SO ₄ (Sulphate)	≤ 10 ppm
As (Arsenic)	≤ 0.1 ppm
Ba (Barium)	≤ 20 ppm
Ca (Calcium)	≤ 10 ppm
Cu (Copper)	≤ 2 ppm
Fe (Iron)	≤ 2 ppm
Mg (Magnesium)	≤ 5 ppm
Na (Sodium)	≤ 0.03 %
Pb (Lead)	≤ 5 ppm



Product Specification

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Jessie Pijck, Head of Laboratory - Haasrode
VWR International bv; Geldenaaksebaan 464; BE-3001
Leuven; Belgium