


AC- 2251xx	CUSTOMER DRAWING	Page 1 of 1
Revision 02	PS Multipurpose Container, 150 ml	 greiner bio-one
	Item-No. 225 1xx	



Cat.-No.	Description	A	B	C	D	E	F
225 170	Multipurpose Beaker	Ø 47.0	Ø 49.0	Ø 54.7	107.0	107.7	10.77
225 180	Multipurpose Beaker	Ø 47.0	Ø 49.0	Ø 54.7	107.0	107.7	10.77

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date May 07, 2007	Date May 07, 2007	Date May 08, 2007	
Date August 1, 2005	Name F. Falter	Name F. Sonntag	Name F. Sonntag	