

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Polygon Spin bar Magnetic Stirring Bars
Country of Origin	Manufactured in United States
Date of Issue (yyyy-mm-dd)	2021-08-26

VWR Catalogue Number	Diameter	Length
74950-284	5 mm ($\frac{3}{16}$ ")	12 mm ($\frac{1}{2}$ ")
74950-286	6 mm ($\frac{1}{4}$ ")	20 mm ($\frac{3}{4}$ ")
74950-288	6 mm ($\frac{1}{4}$ ")	25 mm (1")
74950-290	6 mm ($\frac{1}{4}$ ")	35 mm (1 $\frac{3}{8}$ ")
74950-292	7.9 mm ($\frac{5}{16}$ ")	40 mm (1 $\frac{5}{8}$ ")
74950-294	7.9 mm ($\frac{5}{16}$ ")	44.4 mm (1 $\frac{3}{4}$ ")
74950-296	7.9 mm ($\frac{5}{16}$ ")	51 mm (2")
74950-298	10 mm ($\frac{3}{8}$ ")	60 mm (2 $\frac{5}{16}$ ")
74950-300	10 mm ($\frac{3}{8}$ ")	70 mm (2 $\frac{3}{4}$ ")

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Tolerance +/- .020

Autoclavable @121° C for 20 minutes.

Product Specifications

Material Teflon used to totally encapsulate the magnetic bar meets US Pharmacopeia's Class VI and FDA Regulations 21CFR177.1550, 21CFR177.1520.

Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The product is Latex free.
BPA Statement	Not applicable.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	Not applicable.
Storage Conditions	Store in a clean and dry environment.
Shelf Life	5 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.