

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR® Multiwell cell culture plates
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2025-07-09

VWR Catalogue Number	Type	Description	Culture Area [cm ²]	Recommended Working volume [mL]
732-2905	Sterile	4-well plates, flat bottom, TC-treated	1.96	1.0
732-2904	Sterile	4-well plates, flat bottom, non treated	1.96	1.0
734-2798	Sterile	6-well plates, flat bottom, treated for increased cell attachment	9.6	1.9 – 2.9
734-2323	Sterile	6-well plates, flat bottom, TC-treated	9,6	1.9 – 2.9
734-2777	Sterile	6-well plates, flat bottom, non treated	9.6	1.9 – 2.9
734-2799	Sterile	12-well plates, flat bottom, treated for increased cell attachment	3.85	0.76 – 1.14
734-2324	Sterile	12-well plates, flat bottom, TC-treated	3.85	0.76 – 1.14
734-2778	Sterile	12-well plates, flat bottom, non treated	3.85	0.76 – 1.14
734-2800	Sterile	24-well plates, flat bottom, treated for increased cell attachment	1.93	0.38 – 0.57
734-2325	Sterile	24-well plates, flat bottom, TC-treated	1.93	0.38 – 0.57
734-2779	Sterile	24-well plates, flat bottom, non treated	1.93	0.38 – 0.57
734-2801	Sterile	48-well plates, flat bottom, treated for increased cell attachment	0.83	0.19 – 0.29
734-2326	Sterile	48-well plates, flat bottom, TC-treated	0.83	0.19 – 0.29
734-2780	Sterile	48-well plates, flat bottom, non treated	0.83	0.19 – 0.29
734-2888	Non-Sterile	48-well plates, flat bottom, non treated	0.83	0.19 – 0.29
734-2802	Sterile	96-well plates, flat bottom, treated for increased cell attachment	0.33	0.075 – 0.20

VWR Catalogue Number	Type	Description	Culture Area [cm ²]	Recommended Working volume [mL]
734-2327	Sterile	96-well plates, flat bottom, TC-treated	0.33	0.075 – 0.20
734-2781	Sterile	96-well plates, flat bottom, non treated	0.33	0.075 – 0.20
734-2803	Sterile	96-well plates, round bottom, treated for increased cell attachment	0.58	0.075 – 0.20
734-2328	Sterile	96-well plates, round bottom, TC-treated	0.58	0.075 – 0.20
734-2782	Sterile	96-well plates, round bottom, non treated	0.58	0.075 – 0.20
732-2907	Sterile	384-well plates, flat bottom, TC-treated	0.1135	0.01 – 0.1
732-2906	Sterile	384-well plates, flat bottom, non treated	0.1135	0.01 – 0.1

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 & ISO 13485:2016**

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

The product is Class I category device as defined by FDA in 21CFR Parts 862-892.

Material	Polystyrene (PS).
Sterilization	Product labelled as sterile is EB (Electron beam) irradiated as per ISO 11137 and meets SAL 10 ⁻⁶ with a specified dose range of 15 - 30 kGy which is not applicable for Non-Sterile products.
ATP Assay	Not Applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Testing is conducted to quality all material resins using ISO 10993 standards for cytotoxicity and have been shown to be non-toxic.
Non-Pyrogenic Statement	The acceptance level for product is 0.5 EU/ml or less than 20 EU/device. (TAL Gel Clot Method). Test passed.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not Applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed).
REACH Statement	Not Applicable.
Storage Conditions	Room Temperature.
Shelf Life	Product labelled as sterile has shelf life of 3 years which is not applicable for Non-Sterile products.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	2024-03
2	2.0	Modified Quality system compliance and disclaimer	2025-07
3	1.2		