



Product specification

Product number	2.04706	Version	3.0
Packaging number	2.04706.782	Date of version	05.02.2014
Product name	NBB-Broth-Enrichment-Medium (NBB-B-AM)	Page	1/2

Information about composition and property

Appearance and consistency

reddish, liquid

Odour and taste

unfit for consumption, therefore not applicable

Descriptions

Ingredients

water, beer, sugar, peptone, yeast extract, growth factors, buffer, inhibitors

Application

Microbiological monitoring

Nutrient medium for the detection of beverage and beer spoiling bacteria for the analysis of weak spots within the filling area.

Indication

Please refer to user manual.

Isolation of beer-spoiling bacteria, indicator organism, microbiological monitoring.

pH-value: 5,75 (5,6 - 5,9)

Other information

None.

Storage and packaging information

Packaging and content: Card box, (9 x 250 ml)

Shelf life in the unopened original packaging

At 4-8 °C : 1.080 day(s)

Storage and handling

Do not freeze.

Store under dry and dark conditions.

Chilled storage.

Microbiology

Sterilized and ready for use.

Dangerous good and hazardous material information

No dangerous good.

No hazardous material.



Product specification

Product number	2.04706	Version	3.0
Packaging number	2.04706.782	Date of version	05.02.2014
Product name	NBB-Broth-Enrichment-Medium (NBB-B-AM)	Page	2/2

Disclaimer

This product specification describes the composition and respectively the properties of our product to the best of our knowledge. It replaces all previous versions and is valid without signature.