



Reagecon Red Primary Colour Solution according to European Pharmacopoeia (EP) Chapter 2.2.2

Item No. EPRS01

Reagecon's Red Primary Colour Solution is formulated as a reference solution for the examination of the degree of colouration of liquids in accordance with the European Pharmacopoeia Chapter 2.2.2. The solution may be used as a Primary Colour Solution on its own or as a Standard Colour Solution if used in combination with other Primary Colour Solutions.

Additional Details

Pack Size	100ml
Weight	0.12
Weight Unit	Kg
Colour	Red
Tariff Code	3822 190090
Country of Origin	Ireland

Safety Information

UN Number	1760
Class	8
GHS	

Price €130,⁹⁰ price per 100ml excluding tax

Additional Details

Pack Size	100ml
Weight	0.12
Weight Unit	Kg
Colour	Red
Tariff Code	3822 190090
Country of Origin	Ireland

GHS

DESCRIPTION

This Primary Colour Solution is a ready to use coloration standard presented in high quality, tamper evident bottles for use as a calibration and/or quality control standard. It is prepared with Cobalt Chloride, Hydrochloric Acid and Water.

Manufactured gravimetrically according to European Pharmacopoeia Chapter 2.2.2 guidelines with a Balance Certificate of Calibration issued in accordance with ISO 17025 (INAB ref: 265C)

Concentration verified using a high performance spectrophotometer.

This Primary Colour Solution can be combined in precise ratios to form the Standard Colour Solutions

PRODUCT SPECIFICATIONS

European Pharmacopoeia (EP) Standard, European Pharmacopoeia (EP) Material

Use	Colouration Standard
Tested	Yes
Analyte	Cobalt Chloride

Packaging	Amber Glass Winchester Bottle
Physical Form	Liquid