

Product Specification

Meets ACS Specifications.

Test	Specification
Assay	$\geq 85.5 \%$
Potassium Carbonate (K_2CO_3)	$\leq 0.6 \%$
Chloride (Cl)	$\leq 70 \text{ mg/kg}$
Silver (Ag)	$\leq 0.5 \text{ mg/kg}$
Aluminium (Al)	$\leq 10 \text{ mg/kg}$
Calcium (Ca)	$\leq 10 \text{ mg/kg}$
Iron (Fe)	$\leq 5 \text{ mg/kg}$
Nickel (Ni)	$\leq 3 \text{ mg/kg}$
Pb (Lead)	$\leq 0.4 \text{ mg/kg}$
Sulphate (SO_4)	$\leq 10 \text{ mg/kg}$
Nitrogen Compounds (as N)	$\leq 5 \text{ mg/kg}$
Phosphate (PO_4)	$\leq 4 \text{ mg/kg}$
Silicon Dioxide (SiO_2)	$\leq 5 \text{ mg/kg}$

For Laboratory, Research or Manufacturing use.



Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000