

The VWR bench top transilluminator for stand-alone imaging systems is designed for ease of use.

- Even Illumination- across the filter surface due to special reflector/diffuser to eliminate tube patterns from images- this is essential when images are to be quantified with accuracy.
- High intensity output that equals or betters performance of models with higher wattage tubes.
- Includes a UV protective lid
- Available in one size – 20 x 20cm
- Available with mid wave UV (302nm) excitation wavelength

European Catalogue number	Size	Wavelength
730-1488	20x20cm	302nm