

Axygen® Automation-compatible PCR Microplate

Unique design prevents plate warping and minimizes evaporation during thermal cycling

Warping of PCR microplates through thermal cycling increases sample evaporation and makes them no longer automation-compatible. As a result, Corning invented a unique, one material automation-compatible PCR microplate that sustains flatness post-thermal cycling, has raised-angled rims for improved sealing, and has grid lines for better sample traceability. This microplate is compatible with major automated thermal cyclers, liquid handlers, and sequencers.

Features

- ▶ A rigid frame with a small split **to remain flat** during thermal cycling without distorting its shape.
- ▶ Specially designed raised rims around the wells provide tight sealing to **minimize evaporation**
- ▶ Grid lines, clear wells, and alphanumeric lettering for improved **sample traceability**
- ▶ Compatible with most major thermal cyclers, sequencers, and automated liquid handlers
- ▶ Ultra-thin walls provide fast and consistent heat transfer during PCR.
- ▶ Nonpyrogenic
- ▶ RNase-/DNase-free
- ▶ PCR inhibition and human DNA-free



Ordering Information

Axygen® Automation-compatible PCR Microplate

VWR Cat. No.	Corning Cat. No.	Description	Volume (μL)	Qty/Pk	Qty/Cs
 10015-636	PCR-96-HS-AC-C	Axygen 96-well half skirt automation-compatible, PCR microplate, clear, nonsterile	200	10	50

Accessories

VWR Cat. No.	Corning Cat. No.	Description	Qty/Pk	Qty/Cs
10011-130	AM-96-PCR-RD	Sealing mat for 96-well PCR plate	10	50
47743-943	UC-500	UltraClear pressure sensitive film	100	500
89080-524	PCR-02-FCP-C	0.2 mL clear, thin wall PCR 8-strip flat caps with writing surface, nonsterile	125	1,250
10011-786	PCR-02CP-C	0.2 mL clear, thin wall PCR 8-strip domed caps, nonsterile	125	1,250
47745-512	PCR-2CP-RT-C	0.2 mL ultra-clear, thin wall PCR 8-strip caps for real time PCR, nonsterile	125	1,250
89179-624	PCR-12-FCP-C	0.2 mL clear, thin wall PCR 12-strip flat caps, nonsterile	80	800

Compatible with most major real-time PCR thermal cyclers, including:

- Applied Biosystems 7000, 7300, 7500, 7700, 7900HT, ProFlex™, StepOne™, StepOnePlus™, ViiA™ 7, QuantStudio™ 12K Flex
- BioRad/MJ Research iCycler®, MyiQ™, iQ5™, Chromo4™, DNA Engine Opticon® 2, CFX96 Touch™, CFX Connect™, MiniOpticon™
- Eppendorf Mastercycler® ep realplex

Compatible with most major standard PCR thermal cyclers, including:

- Applied Biosystems Veriti®, 9600, 9700, 2720
- BioRad/MJ Research MyCycler™, T100™, C1000 Touch™, S1000™, DNA Engine Tetrad® 2, PTC-100®, PTC-200®, PTC-225 Tetrad, DNA Engine Dyad®, Dyad Disciple™
- Eppendorf Mastercycler series: ep, pro, nexus, nexus eco, nexus gradient, nexus gradient eco, nexus GSX1, nexus GSX1e, nexus SX1, nexus SX1e

Compatible with the following Applied Biosystems sequencers:

- 3100, 3130, 3700, 3730/3730x

Compatible with most major automated liquid handlers, including:

- Beckman Coulter, Hamilton, Tecan, Agilent, Perkin-Elmer

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

CORNING | FALCON | AXYGEN | GOSSELIN | PYREX



1.800.932.5000 | vwr.com

Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2015 VWR International, LLC. All rights reserved.