

INSTRUCTIONS FOR USE

VWR Catalogue Number	VWRC331942F
Description	Buffer tablets 100 ml pH 7.2
Country of Origin	Manufactured in the UK
Date of Issue (dd/mm/yyyy)	23/11/2020

We herewith state that:

One tablet dissolved in distilled water and made up to 100 ml produces a solution of the stated pH at 20 °C.

The molarity of the buffer solution made up according to these instructions is 0.005 M.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.