

Technical data sheet

Article: Pasteur pipettes 230 mm
Code no.: 926 01 01

Material:
Soda-lime glass

Chemical properties:

- Hydrolytic resistance at 98°C: Class 3 acc. ISO 719
- Hydrolytic resistance at 121°C: Class 3 acc. ISO 720
- Acid resistance: acc. ISO 1776
- Alkali resistance: Class 2 acc. ISO 695

Overall length:
225 - 230 mm

Body length :
85 - 100 mm

Distance of constriction from top :
20 - 30 mm

External diameter tip:
1,2 - 1,8 mm

External diameter body:
6,5 - 8,0 mm

Internal diameter tip:
0,7 -1,5 mm

Hirschmann Laborgeräte
GmbH & Co.KG



Roland Benz
Quality Manager

Eberstadt, 09.01.2023