



Certificate of Conformity

VWR DIGITAL STOPWATCH TRACEABLE®

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat No: 609-0212

Range: Timing range is 9 hours, 59 minutes, 59 seconds, and 99 hundredths

Resolution: 1/100 second

Colour: Black/blue

Size (WXDXH): 63.5 x 19 x 76.2 mm

Other Characteristics:

- Quartz-crystal stopwatch with large LCD digits.
- Functions: Single action, time out, cumulative split, interval split, continuous roll-over timing, countdown alarm timing.
- Tough, water-resistant ABS plastic casing seals out moisture, dust and fumes.
- Ideal for enzyme tests, viscosity measurements, animal studies, and all timed lab tests.
- In Compliance with ISO/IEC 17025; accredited by A2LA indicating traceability to standards provided by NIST.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **April 2018**