



Certificate of Conformity

VWR Mat Frame for Purestep Mat

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

VWR Catalog No.	L×W	For mats L×W	Colour	Thickness	Pk
115-1287	508×965 mm	457×914 mm	White	3.2 mm	1
115-1288	508×1194 mm	457×1143 mm	White	3.2 mm	1
115-1289	660×965 mm	609×914 mm	White	3.2 mm	1
115-1290	711×1193 mm	660×1143 mm	White	3.2 mm	1
115-1291	965×965 mm	914×914 mm	White	3.2 mm	1
115-1292	965×1574 mm	914×1524 mm	White	3.2 mm	1
115-1296	965×1194 mm	914×1140 mm	White	3.2 mm	1

Material: Polystyrene

Other Characteristics:

- Protect floors from adhesive residue.
- Reduces the risk of slips and falls.
- Non slip parts made from polymer coated, polyester textile.
- Naturally antimicrobial, mould and mildew resistant.

This document has been produced electronically and is valid without a signature

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **August 2019**