



Specification

1.01299.0001 Enterococci 100 Readycult®

Specification

Appearance (clearness)	clear
Appearance (color)	yellowish to yellowish-brown
pH-value (25 °C)	7.3 - 7.7
Sterility test (48 hrs., 35 °C)	passes test

Specification

Inoculum on reference medium (Enterococcus faecalis ATCC 19433 (WDCM 00009))	10 - 100
Inoculum on reference medium (Enterococcus durans ATCC 19432)	10 - 100
Inoculum on reference medium (Staphylococcus aureus ATCC 25923 (WDCM 00034))	≥ 1000
Growth (Enterococcus faecalis ATCC 19433 (WDCM 00009))	fair to good
Growth (Enterococcus durans ATCC 19432)	fair to good
Growth (Staphylococcus aureus ATCC 25923 (WDCM 00034))	fair to good
Color change to green-blue (Enterococcus faecalis ATCC 19433 (WDCM 00009))	+
Color change to green-blue (Enterococcus durans ATCC 19432)	+
Color change to green-blue (Staphylococcus aureus ATCC 25923 (WDCM 00034))	-

Incubation: 24 - 48 hrs.; 35 °C; aerobic.

Stefanie Fischer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.