



Product Specification Sheet

VWR CRITICAL SWAB, DOUBLE, HIGHLY ABSORBENT COTTON HEADS, PACKAGED IN AN ANTI-STATIC LEAF PACK, WHITE

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Cat. N°: **149-0341**

Dimensions:

- **Handle length (mm) : 78.0**
- **Handle width (mm) : 2.5**
- **Head length (mm) : 15.0**
- **Head width (mm) : 4.4**
- **Total length (mm) : 93.0**

Qualities of swab :

- **Easy disposal**
- **Lint-less and residue free**
- **Highly absorbent cotton does not swell in solvents**
- **Low in Non-volatile residues**
- **Burnable**
- **Environmental friendly**



VWR International bvba/sprl, Haasrode Research Park Zone 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>

Version: 2.0



Product Specification Sheet

VWR CRITICAL SWAB, DOUBLE, HIGHLY ABSORBENT COTTON HEADS, PACKAGED IN AN ANTI-STATIC LEAF PACK, WHITE

Applications:

- Health and beauty
- Cleaning disk drives, optics, PC board assembly and other critical industrial components
- Cleaning sensitive components
- Veterinary
- Weapons cleaning

Non-Volatile Residue (NVR):

NVR

DI	0.155 mg/tip
IPA	0.102 mg/tip
Acetone	0.112 mg/tip
Hexane	Not detected

Ion Chromatography Testing Results:

Ionic Cleanliness

Anions

Chloride 1.4.39 µg / tip

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there's any doubt ask us.

Date of issue: **September 2012**

