

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Flat Caps for Screw Cap Microtubes
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2025-02-26

VWR Catalogue Number	Cap Colour	Sterility	Package
76417-216	Clear	Non-sterile	Pack of 500
			Case of 5000
76417-230	Opaque assortment	Non-sterile	Pack of 500
			Case of 5000
76417-218	Opaque blue	Non-sterile	Pack of 500
			Case of 5000
76417-220	Opaque green	Non-sterile	Pack of 500
			Case of 5000
76417-222	Opaque orange	Non-sterile	Pack of 500
			Case of 5000
76417-224	Opaque red	Non-sterile	Pack of 500
			Case of 5000
76417-226	Opaque violet	Non-sterile	Pack of 500
			Case of 5000
76417-228	Opaque yellow	Non-sterile	Pack of 500
			Case of 5000
76437-218	White	Non-sterile	Pack of 500
			Case of 5000

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Polypropylene (PP); Compliance with FDA 21 CFR 177.1520 & USP Class VI standards.
Sterilization	The products are non-sterile.
ATP Assay	The Products are free from ATP Assay.
DNase & RNase Free	The products are free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Testing is conducted to quality all material resins and meet the requirements of a USP Class VI standards and have been shown to be non-toxic.
Non-Pyrogenic Statement	The acceptance level for product is less than 0.06 EU/ml. (Limulus Amoebocyte Lysate (LAL) assay protocol and USP guidelines.
Latex Statement	The products are Latex free.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr ⁶⁺), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.

REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Store at Normal room temperature.
Shelf-Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.

This document has been produced electronically and is valid without a signature.