

Certificate of Conformity

VWR Catalogue Number	215-4497
Description	HDPE Bottle with Plug, 150 mL capacity, square bottle, sterile, red cap, 45/400 mm neck size and finish, 50.5 mm length, 50.5 mm width
Lot Number	2510090865
Date of Sterilization (yyyy-mm-dd)	2025-10
Expiry Date (yyyy-mm-dd)	2027-09
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2026-06-08

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 & ISO 15378:2017**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Test Name	Test Reference	Accuracy	Result
Leak Test	SOP/QA/07A	General Insp. Level-1, AQL -1.5	Leak Proof
Ethylene Oxide Dose	SOP/PRO/21	NA	Pass
Dimensional	SOP/QA/01	General Insp. Level-1, AQL -1.5	Within Tolerance

Product Specifications

Material

Bottle: High Density Polyethylene, 100% Virgin, FDA 21 CFR Compliance, USP Class VI Passed, Colour: Natural.

Closure with Plug: Polypropylene, 100% Virgin, FDA 21 CFR Compliance, USP Class VI Passed.

Monomers and additives used to produce the products are listed in List of Authorised Substances of Regulation EU/10/2011.

VWR International bv, Haasrode Researchpark 2020, Geldenaaksebaan 464, B-3001 Leuven, Belgium - <http://www.vwr.com>

The material grades conform to EU Directive 1935/2004 on materials and articles intended to come into contact with food.

Sterilization	Product labelled as sterile is ETO (Ethylene Oxide) irradiated as per ISO 11135:2014 and meets SAL 10 ⁻³ .
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	This product is free of cytotoxicity effects.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Nitrosamines	The raw material used to manufacture above product does not contain Nitrosamines.

Melamine	Melamine is not intentionally used and is not expected to be present in the product based on the composition of the raw materials used in the manufacture of the product.
Storage Conditions	Store at 25 + 3 / - 8 °C, Relative Humidity 50 ± 10 % RH.
Shelf Life	2 years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.