

VWR Brand Autoclave Bags

VWR CATALOG NUMBERS:

[Manufacturing country of origin for all below: USA]

11215-822	11215-823	11215-824	11215-826	14232-184	14232-764
-----------	-----------	-----------	-----------	-----------	-----------

DESCRIPTION:

Double Thick biohazard bags made from autoclavable High Density Polyethylene film, specifically engineered to exhibit seal ability, good toughness and autoclavable temperature resistance (typical program of 20 mins. @ Operating temperature to 250°F, film will melt above 260°F) as utilized in some packaging and sterilization applications.

Meets Dart Impact Strength requirements for all 50 States in USA.

TYPICAL PRODUCT CHARACTERISTICS:

Values expressed below for available nominal film thicknesses.

PROPERTIES		4.0 mil	UNITS	TEST METHOD
Elmendorf Tear Resistance	MD TD	260 325	grams grams	ASTM D1922
Dart Impact Strength		400	grams	ASTM D1709
Coefficient of Friction		<0.30		ASTM D1894
Yield		7284	sq.in./lb.	Tufpak

FOR YOUR PROTECTION:

The information and recommendations in this publication are, to the best of our knowledge, accurate and reliable. Suggestions made concerning uses or applications are only the opinion of Tufpak, Inc., and users should make their own tests to determine the suitability of these products for their own particular purposes. However, because of numerous factors affecting results, Tufpak, Inc. makes no warranty of any kind, express or implied, including those of merchantability and fitness for purpose, other than that the material conforms to its Standard Specifications. Statements herein, therefore, should not be construed as representations or warranties. The responsibility of Tufpak, Inc. for claims arising out of breach of warranty, negligence, strict liability, or otherwise is limited to the purchase price of the material. Statements concerning the use of products described herein are not to be construed as recommending the infringement of any patent and no liability for infringement arising out of any such use is assumed. 1/2021