

Orbi-Shaker Jr.	
470318-958	<i>Orbi-Shaker™Jr</i> with rubber mat platform (9.5"x8"), 100-240V (US Plug)
470318-960	<i>Orbi-Shaker™Jr</i> with rubber mat platform (9.5"x8"), 100-240V (European 2 prong plug)
Accessories	
470200-848	<i>MAGic Clamp™</i> universal platform for flasks & tube racks (11x9.5in)
470318-982	<i>MAGic Clamp™</i> magnetic clamp, 1000ml Erlenmeyer (max. 4)
470200-850	<i>MAGic Clamp™</i> magnetic clamp, 500ml Erlenmeyer (max. 6)
470200-852	<i>MAGic Clamp™</i> magnetic clamp, 250ml Erlenmeyer (max. 8)
470200-854	<i>MAGic Clamp™</i> magnetic clamp, 125ml Erlenmeyer (max. 12)
470017-186	<i>MAGic Clamp™</i> magnetic clamp, 50ml Erlenmeyer (max. 20)
470318-986	<i>MAGic Clamp™</i> magnetic clamp, 25ml Erlenmeyer
470318-992	<i>MAGic Clamp™</i> magnetic clamp, one microplate (max. 4)
470200-856	<i>MAGic Clamp™</i> Tube Rack, 30x15ml & 20x50ml tubes (max. 1)
470318-996	<i>MAGic Clamp™</i> Tube Rack, 32x15ml, tubes (max.2)
470319-002	<i>MAGic Clamp™</i> Tube Rack, 15x50ml, tubes (max. 2)
470318-990	<i>MAGic Clamp™</i> COMBO, for 125ml, 250ml, 500ml flasks
470319-008	<i>MAGic Clamp™</i> tilted holder for 15mL and 50mL tube racks (max 2)

IMPORTANT NOTE

Orbi-Shakers come with a standard platform with rubber mat. If you want to use MAGic Clamp™ system, you must also purchase the MAGic Clamp™ Universal Platform recommended for your model. MAGic Clamp™. The standard rubber mat platforms are not magnetic and wil not work with the MAGic Clamp™ System.

Orbi-Shaker - Orbi Shaker XL - Orbi-Shaker CO2	
470318-948	<i>Orbi-Shaker™XL</i> orbital shaker with flat mat platform, 115V
470318-950	<i>Orbi-Shaker™XL</i> orbital shaker with flat mat platform, 230V
470318-954	<i>Orbi-Shaker™</i> with rubber mat platform (13"x12") 115V
470318-956	<i>Orbi-Shaker™</i> with rubber mat platform (13"x12"), 230V
470318-962	<i>Orbi-Shaker™ CO2 with remote controller & rubber mat platform (13"x12"), 115V</i>
470318-964	<i>Orbi-Shaker™ CO2 with remote controller & rubber mat platform (13"x12"), 230V</i>
Accessories	
470200-848	<i>MAGic Clamp™</i> universal platform (SM) for flasks & tube racks (9.5 x 11.5")
470318-952	<i>MAGic Clamp™</i> universal platform (LG) for flasks & tube racks, 14 x 12"
470318-982	<i>MAGic Clamp™</i> magnetic clamp, 1000ml Erlenmeyer (max. 5)
470200-850	<i>MAGic Clamp™</i> magnetic clamp, 500ml Erlenmeyer (max. 9)
470200-852	<i>MAGic Clamp™</i> magnetic clamp, 250ml Erlenmeyer (max. 14)
470200-854	<i>MAGic Clamp™</i> magnetic clamp, 125ml Erlenmeyer (max. 24)
470017-186	<i>MAGic Clamp™</i> magnetic clamp, 50ml Erlenmeyer (max. 30)
470318-986	<i>MAGic Clamp™</i> magnetic clamp, 25ml Erlenmeyer
470318-992	<i>MAGic Clamp™</i> magnetic clamp, one microplate (max. 6)
470200-856	<i>MAGic Clamp™ Tube Rack, 30x15ml & 20x50ml tubes (max. 2)</i>
470318-994	<i>MAGic Clamp™ Tube Rack, 72x14mm culture tubes</i>
470318-996	<i>MAGic Clamp™</i> Tube Rack, 32x15ml, tubes (max. 3)
470319-002	<i>MAGic Clamp™</i> Tube Rack, 15x50ml, tubes (max. 3)
470319-006	<i>MAGic Clamp™</i> Tube Rack, 4x500ml or 600ml conical bottles (max. 1)
470319-000	<i>MAGic Clamp™</i> Tube Rack, 12 x 1.5/2.0ml tubes, horizontal
470318-998	<i>MAGic Clamp™</i> Tube Rack, 4 x 15ml tubes, horizontal
470319-004	<i>MAGic Clamp™</i> Tube Rack, 2 x 50ml tubes, horizontal
470319-008	<i>MAGic Clamp™</i> tilted holder for 15mL and 50mL tube racks (max 3)
470318-990	<i>MAGic Clamp™</i> COMBO, for 125ml, 250ml, 500ml flasks
470319-014	Dedicated Platform, 5 x 500ml Erlenmeyer flasks
470319-012	Dedicated Platform, 8 x 250ml Erlenmeyer flasks
470319-010	Dedicated Platform, 12 x 125ml Erlenmeyer flasks
470319-016	Dedicated Platform, holds 4 standard micro plates
470319-018	Universal Spring Platform