

# Certificate of Conformity

|                                   |                                                                                                                                                  |
|-----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>VWR Catalogue Number</b>       | 548-0048A                                                                                                                                        |
| <b>Description</b>                | Shell vials, Vials, Flat bottom base, Borosilicate glass 5.1, No cap included, 15 mm PE plug, 4 mL Capacity, 15 mm diameter, 45 mm height, Clear |
| <b>Lot Number</b>                 | 000000158309                                                                                                                                     |
| <b>Country of Origin</b>          | Manufactured in Germany                                                                                                                          |
| <b>Date of Issue (yyyy-mm-dd)</b> | 2023-06-23                                                                                                                                       |

## Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

## Product Specifications

|                                |                                                                                                                                                                                                                                                                                                                                                              |
|--------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Material</b>                | Borosilicate glass 5.1.                                                                                                                                                                                                                                                                                                                                      |
| <b>Sterilization</b>           | Not applicable.                                                                                                                                                                                                                                                                                                                                              |
| <b>ATP Assay</b>               | Not applicable.                                                                                                                                                                                                                                                                                                                                              |
| <b>DNase &amp; RNase Free</b>  | Not applicable.                                                                                                                                                                                                                                                                                                                                              |
| <b>BSE/TSE</b>                 | No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE). |
| <b>Cytotoxicity</b>            | Not applicable.                                                                                                                                                                                                                                                                                                                                              |
| <b>Non-Pyrogenic Statement</b> | Not applicable.                                                                                                                                                                                                                                                                                                                                              |

|                           |                                                                                                                                                                                                                 |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Latex Statement</b>    | Not applicable.                                                                                                                                                                                                 |
| <b>BPA Statement</b>      | Not applicable.                                                                                                                                                                                                 |
| <b>DEHP Statement</b>     | Not applicable.                                                                                                                                                                                                 |
| <b>RoHS</b>               | Not applicable.                                                                                                                                                                                                 |
| <b>REACH Statement</b>    | No substances as listed on >Substance of Very High Concern< candidate list published 2023-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed. |
| <b>Storage Conditions</b> | Not applicable.                                                                                                                                                                                                 |
| <b>Shelf Life</b>         | Not applicable.                                                                                                                                                                                                 |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature

## Version History

| Sl.no | Version number | Description                          | Month & Year |
|-------|----------------|--------------------------------------|--------------|
| 1     | 1.0            | Created General CoQ using new layout | 12-2022      |
|       |                |                                      |              |
|       |                |                                      |              |