

PRODUCT DATA SHEET

Combi-Stopper

Universal Closing Stopper



Combi-Stoppers are closing cones with luer lock fitting male and female. They are universally suitable for sterile closing of all open male or female luer lock and luer slip connections and provide a contact-free luer taper.

- The Combi-Stoppers are made of Polyethylene
- For closing infusion systems, e.g., syringes, IV-sets or IV-catheters
- Can be used for male and female luer accesses.

Advantages

Easy access for Luer-Slip and Luer-Lock connections

Sterile closing of all kinds of open IV-accesses in IV-Sets and of pre-filled syringes.

Material

- Not manufactured with BPA
- Not manufactured with DEHP
- Not manufactured with Latex / Natural rubber
- Not manufactured with PVC

Compatibility

- Combi-Stoppers can be used with alcohol and lipids
- Combi-Stoppers are non-toxic and non-pyrogenic.

Pressure

Combi-Stoppers are pressure resistant up to 2 bar.

Portfolio

- Available in red, blue and white
- Box of 100 pcs

Pediatric use

The device can be used for all patients for which infusion therapy is prescribed. No gender or age-related limitations. Combi-Stoppers can be used for adults, pediatric and neonates.

Article Description	Article Number
COMBI-STOPPER RED	4495101
COMBI-STOPPER RED	4495101R
COMBI-STOPPER BLUE	4495152
COMBI-STOPPER BLUE	4495152B
COMBI-STOPPER WHITE	4495209