



Material	83872.270
Material description	Nitric acid 67-69%
Grade	NORMATOM for trace metal analysis
CAS Number	7697-37-2
Molecular formula	HNO ₃
Molecular mass	63.01

Characteristics	Specifications
Assay	67 - 69 %
Colouration	≤ 10 APHA
Cl (Chloride)	≤ 0.2 ppm
Total P (Phosphorus)	≤ 0.01 ppm
Total S (Sulphur)	≤ 0.3 ppm
Ag (Silver)	≤ 0.1 ppb
Al (Aluminium)	≤ 1 ppb
As (Arsenic)	≤ 0.5 ppb
Au (Gold)	≤ 0.1 ppb
B (Boron)	≤ 1 ppb
Ba (Barium)	≤ 0.1 ppb
Be (Beryllium)	≤ 0.1 ppb
Bi (Bismuth)	≤ 0.1 ppb
Ca (Calcium)	≤ 1 ppb
Cd (Cadmium)	≤ 0.5 ppb
Ce (Cerium)	≤ 0.1 ppb
Co (Cobalt)	≤ 0.5 ppb
Cr (Chromium)	≤ 1 ppb
Cs (Cesium)	≤ 0.1 ppb
Cu (Copper)	≤ 0.5 ppb
Dy (Dysprosium)	≤ 0.1 ppb
Er (Erbium)	≤ 0.1 ppb
Eu (Europium)	≤ 0.1 ppb
Fe (Iron)	≤ 1 ppb
Ga (Gallium)	≤ 0.1 ppb
Gd (Gadolinium)	≤ 0.1 ppb
Ge (Germanium)	≤ 0.1 ppb
Hf (Hafnium)	≤ 0.1 ppb

>>> Continued on page 2 >>>

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.



Characteristics	Specifications
Hg (Mercury)	≤ 0.1 ppb
Ho (Holmium)	≤ 0.1 ppb
In (Indium)	≤ 0.1 ppb
K (Potassium)	≤ 1 ppb
La (Lanthanum)	≤ 0.1 ppb
Li (Lithium)	≤ 0.1 ppb
Lu (Lutetium)	≤ 0.1 ppb
Mg (Magnesium)	≤ 1 ppb
Mn (Manganese)	≤ 0.1 ppb
Mo (Molybdenum)	≤ 0.1 ppb
Na (Sodium)	≤ 1 ppb
Nb (Niobium)	≤ 0.1 ppb
Nd (Neodymium)	≤ 0.1 ppb
Ni (Nickel)	≤ 0.5 ppb
Pb (Lead)	≤ 0.1 ppb
Pd (Palladium)	≤ 0.5 ppb
Pr (Praseodymium)	≤ 0.1 ppb
Pt (Platinum)	≤ 0.5 ppb
Rb (Rubidium)	≤ 0.1 ppb
Re (Rhenium)	≤ 0.1 ppb
Rh (Rhodium)	≤ 0.5 ppb
Ru (Ruthenium)	≤ 0.5 ppb
Sb (Antimony)	≤ 0.5 ppb
Sc (Scandium)	≤ 0.1 ppb
Se (Selenium)	≤ 1 ppb
Sm (Samarium)	≤ 0.1 ppb
Sn (Tin)	≤ 0.5 ppb
Sr (Strontium)	≤ 0.1 ppb
Ta (Tantalum)	Information only
Tb (Terbium)	≤ 0.1 ppb
Te (Tellurium)	≤ 0.1 ppb
Th (Thorium)	≤ 0.1 ppb
Ti (Titanium)	≤ 0.5 ppb
Tl (Thallium)	≤ 0.1 ppb
Tm (Thulium)	≤ 0.1 ppb
U (Uranium)	≤ 0.1 ppb
V (Vanadium)	≤ 0.5 ppb
W (Tungsten)	≤ 0.1 ppb
Y (Yttrium)	≤ 0.1 ppb
Yb (Ytterbium)	≤ 0.1 ppb
Zn (Zinc)	≤ 0.5 ppb
Zr (Zirconium)	≤ 0.1 ppb

>>> Continued on page 3 >>>



Product Specification



Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International

