

Combination of epT.I.P.S.® Pipette Tips with Research® plus Pipettes

Eppendorf Research® plus	0.1 – 10 µL dark gray	0.1 – 10 µL M medium gray	0.1 – 20 µL medium gray	0.5 – 20 µL L light gray	0.1 – 20 µL pearl white	0.5 – 100 µL light yellow	2 – 20 µL yellow	2 – 100 µL yellow	2 – 200 µL yellow	20 – 300 µL orange	50 – 1,000 µL blue	50 – 1,250 µL green	50 – 1,250 µL L dark green	0.25 – 2.5 mL red	0.1 – 5 mL violet	0.2 – 5 mL L violet	0.5 – 10 mL turquoise	0.5 – 10 mL L turquoise	GEloader™/ Microloader™ 0.5-20 µL
epT.I.P.S.®	■	■	■	■	□	■	■	■	■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S.® LoRetention	■			■					■		■		■				■	■	■
ep Dualfilter T.I.P.S.®	■	■		■	■	■	■	■	■	■	■		■		■	■	■	■	
ep Dualfilter T.I.P.S.® LoRetention	■			■				■		■	■								
ep Dualfilter T.I.P.S.® SealMax	■			■				■	■	■	■								
Fixed volume, single channel pipettes																			
10 µL	■	■	■	■															■
20 µL	■	■	■	■															■
10 µL	■						■	■	■										
20 µL	■						■	■	■	■									
25 µL/50 µL/100 µL	■							■	■	■									
200 µL	■								■	■									
200/250/500/1,000 µL	■										■	■	■						
Variable volume, single-, 8- and 12-channel pipettes																			
0.1 – 2.5 µL	■	■	■	■															
0.5 – 10 µL	■	■	■	■															■
2 – 20 µL	■	◆	◆	■	■														■
2 – 20 µL	■						■	■	■	◆									
10 – 100 µL	■						◆	■	■	■									
20 – 200 µL	■						◆	◆	■	■									
30 – 300 µL	■						◆	◆	◆	■									
100 – 1,000 µL	■										■	■	■						
120 – 1,200 µL	■											■	■						
0.25 – 2.5 mL	■													■		■			
0.5 – 5 mL	■														■	■			
1 – 10 mL	■																■	■	
Variable volume, 16- and 24-channel pipettes																			
1 – 20 µL	□				■														
5 – 100 µL	■					■													
Move It™, variable volume, 4, 6, 8, and 12-channel pipettes																			
1 – 20 µL	□				■														
5 – 100 µL	■					■													
15 – 300 µL	■						◆	◆	◆	■									
50 – 1,200 µL	■											■	■						
■ normal volume ◆ limited volume																			

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf®, the Eppendorf Brand Design, epT.I.P.S.®, ep Dualfilter T.I.P.S.® and Eppendorf Research® are registered trademarks of Eppendorf AG, Germany.
GEloader™, Microloader™ and Move It™ are trademarks of Eppendorf AG, Germany.
All rights reserved, including graphics and images. Copyright © 2019 by Eppendorf AG.

For more information go to
www.eppendorf.com/eptips

Combination of epT.I.P.S.® Pipette Tips with Reference® 2 Pipettes

Eppendorf Reference® 2	0.1 – 10 µL S dark gray	0.1 – 10 µL M medium gray	0.1 – 20 µL medium gray	0.5 – 20 µL L light gray	2 – 20 µL yellow	2 – 100 µL yellow	2 – 200 µL yellow	20 – 300 µL orange	50 – 1,000 µL blue	50 – 1,250 µL green	50 – 1,250 µL L dark green	0.25 – 2.5 mL red	0.1 – 5 mL violet	0.2 – 5 mL L violet	0.5 – 10 mL turquoise	0.5 – 10 mL L turquoise	GEloader™/ Microloader™ 0.5-20 µL
	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S.®	■		■	■			■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S.® LoRetention	■			■			■		■				■	■		■	
ep Dualfilter T.I.P.S.®	■	■		■	■	■	■	■	■		■		■	■	■		
ep Dualfilter T.I.P.S.® LoRetention	■			■		■		■	■								
ep Dualfilter T.I.P.S.® SealMax	■			■		■	■	■	■								
Fixed volume																	
1 µL/2 µL	■	■	■														■
5 µL/10 µL	■	■	■	■													■
20 µL	■		■	■													
10 µL	■				■	■	■										
20 µL	■				■	■	■	■									
25 µL/50 µL/100 µL	■					■	■	■									
200 µL	■						■	■									
200/250/500/1,000 µL	■								■	■	■						
2 mL/2.5 mL	■											■					
Adjustable volume																	
0.1 – 2.5 µL	■	■	■	■													
0.5 – 10 µL	■	■	■	■													■
2 – 20 µL	■	◆	◆	■													■
2 – 20 µL	■				■	■	■	◆									
10 – 100 µL	■				◆	■	■	■									
20 – 200 µL	■				◆	◆	■	■									
30 – 300 µL	■				◆	◆	◆	■									
100 – 1,000 µL	■								■	■	■						
0.25 – 2.5 mL	■											■					
0.5 – 5 mL	■												■	■			
1 – 10 mL	■														■	■	

■ normal volume ◆ limited volume

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf®, the Eppendorf Brand Design, epT.I.P.S.®, ep Dualfilter T.I.P.S.® and Eppendorf Reference® are registered trademarks of Eppendorf AG, Germany.
GEloader™ and Microloader™ are trademarks of Eppendorf AG, Germany.
All rights reserved, including graphics and images. Copyright © 2019 by Eppendorf AG.

For more information go to
www.eppendorf.com/eptips

Combination of epT.I.P.S.® Pipette Tips with Xplorer®/Xplorer® plus Pipettes

Eppendorf Xplorer®/Eppendorf Xplorer® plus	0.1 – 10 µL dark gray	0.1 – 10 µL M medium gray	0.1 – 20 µL medium gray	0.5 – 20 µL L light gray	0.1 – 20 µL pearl white	0.5 – 100 µL light yellow	2 – 20 µL yellow	2 – 100 µL yellow	2 – 200 µL yellow	20 – 300 µL orange	50 – 1,000 µL blue	50 – 1,250 µL green	50 – 1,250 µL L dark green	0.25 – 2.5 mL red	0.1 – 5 mL violet	0.2 – 5 mL L violet	0.5 – 10 mL turquoise	0.5 – 10 mL L turquoise	GEloader™/ Microloader™ 0.5 – 20 µL
	■	■	■	■	□	■	■	■	■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S.®	■		■	■	■	■			■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S.® LoRetention	■			■					■		■				■	■			
ep Dualfilter T.I.P.S.®	■	■		■	■	■	■	■	■	■	■		■		■	■	■	■	
ep Dualfilter T.I.P.S.® LoRetention	■			■				■		■	■								
ep Dualfilter T.I.P.S.® SealMax	■			■				■	■	■	■								
Variable volume, single-, 8- and 12-channel pipettes																			
0.5 – 10 µL	■	■	■	■															■
1 – 20 µL	■	■	■	■															
5 – 100 µL	■						■	■	■	■									
10 – 200 µL	■						■	■	■	■									
15 – 300 µL	■						■	■	■	■									
50 – 1,000 µL	■										■	■	■						
50 – 1,200 µL	■											■	■						
0.125 – 2.5 mL	■													■					
0.2 – 5 mL	■														■	■			
0.5 – 10 mL	■																■	■	
Variable volume, 16- and 24-channel pipettes																			
1 – 20 µL	□				■														
5 – 100 µL	■					■													
Move It™, variable volume, 4, 6, 8, and 12-channel pipettes																			
1 – 20 µL	□				■														
5 – 100 µL	■					■													
15 – 300 µL	■						■	■	■	■									
50 – 1,200 µL	■											■	■						
■ normal volume ♦ limited volume																			

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

Eppendorf®, the Eppendorf Brand Design, epT.I.P.S.®, ep Dualfilter T.I.P.S.® and Eppendorf Xplorer® are registered trademarks of Eppendorf AG, Germany.
GEloader™, Microloader™ and Move It™ are trademarks of Eppendorf AG, Germany.
All rights reserved, including graphics and images. Copyright © 2019 by Eppendorf AG.

For more information go to
www.eppendorf.com/eptips