



Product Specification Sheet

VWR REPLACEMENT SPONGES FOR CLEANING SYSTEMS (A)

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Description	Pk
129-2506	General purpose sponge with stainless steel clip (254×83 mm), blue	12
129-2507	General purpose sponge with stainless steel clip (355×83 mm), blue	12
129-2505	General purpose sponge with galvanised clip (355×83 mm), blue	12
129-2508	White polyester sponge with stainless steel clip (355×83 mm)	12
129-2509	Microfibre sponge with stainless steel clip (355×83 mm)	12

Material: Sponge: Hydrophilic Ether Foam, covered or not with polyester or microfibre (according to model)
Channel: Stainless steel or Galvanised steel (according to model)

Other Characteristics:

- These systems can be used without a wringer.
- They are particularly suitable for applying disinfectants to floors, walls, and ceilings.
- The new, high performance foam for clean rooms enables excellent and lasting liquid transfer.
- Autoclavable
- Resistant to gamma rays and ETO(Ethylene Oxide Sterilization)
- Resistant to most widely available disinfectants and detergents
- Sponges made from laminated polyester for liquid absorption and highly delicate, lint-free disinfection
- Sponges made from laminated microfibres for liquid absorption, deep cleaning and highly delicate, lint-free disinfection

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **January 2014**