




Time® Tape The Toughest, Most Reliable Labeling Tape in the Industry!

All purpose, self-adhesive Time® Tape is well-suited for a wide variety of labeling applications. It can be written on with pen, pencil, or permanent marker and has a removable adhesive. This strong, coated paper tape is oil, water, and acid resistant. It will withstand temperatures from -23°C to +121°C and has an application temperature as low as 0°C.

Stock Time® Tape 500" (12,7 m) Rolls

Stock Time® Tape in 12,7 m Rolls - Available in the 18 stock colours shown below			
Product Code	Tape Width	Core Size	Qty.
T-512	1,3 cm	2,5 cm	1270 cm/RL
T-534	1,9 cm		
T-501	2,5 cm		
T-5112	3,8 cm		
T-502	5 cm		

Stock Time® Tape 2160" (54,86 m) Rolls

Stock Time® Tape in 54,86 m Rolls - Available in the 18 stock colours shown below			
Product Code	Tape Width	Core Size	Qty.
T-1260	1,3 cm	7,6 cm	5486 cm/RL
T-3460	1,9 cm		
T-160	2,5 cm		
T-11260	3,8 cm		
T-260	5 cm		

Time® Tape, 500" (12,7 m) Rainbow Packs

Time® Tape is also available in a convenient Rainbow pack, which contains assorted rolls of 500" (12,7 m) tape. The Rainbow pack is available in the three widths below and priced per box. Order in full box quantities.

Stock Time® Tape in 500" (12,7 m) Rolls - Rainbow Pack			
Product Number	Tape Width	Colours per Box	Box Qty.
T-512-RP	1,3 cm	(5) White, (4) Yellow, (3) Blue, (3) Green, (3) Orange, (3) Red, (3) Pink	24 Rolls
T-534-RP	1,9 cm	(4) White, (2) Yellow, (2) Blue, (2) Green, (2) Orange, (2) Red, (2) Pink	16 Rolls
T-501-RP	2,5 cm	(2) White, (2) Yellow, (2) Blue, (2) Green, (1) Orange, (2) Red, (1) Pink	12 Rolls

COLOUR CHART

1 - White	8 - Tan	14 - Grey
2 - Yellow	9 - Salmon	15 - Lime
3 - Green	10 - Copper	16 - Aqua
4 - Pink	11 - Rose	17 - Chartreuse
5 - Red	12 - Violet	18 - Black
6 - Orange	13 - Lavender	19 - Fl. Red (Custom)
7 - Blue		



Custom Time® Tape Options



When ordering custom imprinted Time® Tape, please specify the following:

1) Length of the roll and quantity you need:

• 12,7 m (24 roll min) • 54,86 m (6 roll min)

2) Width of the roll you need:

1,3 cm, 1,9 cm, 2,5 cm, 3,8 cm, or 5 cm

3) Repeat length of the imprint:

Repeat length is the distance from the beginning of one imprint to the beginning of the next imprint.

REPEAT REPEAT REPEAT

(The repeat length = 3,8 cm)

4) Exact wording needed:

For best results, include a sample or a sketch.

5) Tape Colour:

Choose from one of the 19 colours (shown on page 6).

6) Colour Matching:

Any colour needed, other than the 19 shown is subject to a \$25.00 ink matching charge. Add 20% for each additional imprint ink colour.

7) Reverse Plates:

When a light coloured type on a dark background is needed, there is an additional 25% tinting charge.

8) Perforations:

Add an additional 10% to the tape price for perforations between imprints.


9) Special Logos:

Custom Logos or unique graphics/type are subject to an additional plate charge of \$35.00 for one colour and \$50.00 for two colours.

10) Proof:

Upon written request with your purchase order, we will supply a one colour proof of your label. This proof must be returned within 30 days with your approval before the order can be processed. Proof orders that are cancelled will be billed for all production charges incurred prior to cancellation.

Single and Multiple Roll Plastic Dispensers




HD-2



1052P-B




500P




500P-LBL-INS



60P



3161



500P-LAN



500P-LAN-LBL

Product Number	Description
HD-2	Hand-held tape and label dispenser. Holds a variety of tape and/or label products up to a total of 2,8 cm. Includes 2 spacers and a plastic clip to lock dispenser. Can be used as a single or multiple roll dispenser with pre-dispensing label mode.
1052P-B	Can hold multiple rolls of MedVision® products up to a total of 21,6 cm. Includes inserts, hanging bars and bottom suction cups for secure placement. Size: (H) 11,7 cm x (W) 22,8 cm x (D) 15 cm
500P	For rolls with 2,5 cm cores. Can hold multiple rolls of tape and label products up to a total of 22,8 cm. Includes flip-top cover and bottom suction cups for secure placement. Size: (H) 8,6 cm x (W) 25 cm x (D) 11,4 cm
500P-LBL-INS	Insert slots for the 500P dispenser. Holds up to (11) 1,3 cm rolls.
60P	For rolls with 7,6 cm cores. Can hold multiple rolls of tape and label products up to a total of 23,5 cm. Includes flip-top cover and bottom suction cups for secure placement. Size: (H) 15,2cm x (W) 25 cm x (D) 18,4 cm
3161	Weighted dispenser. Adjustable for rolls with 2,5 cm or 7,6 cm cores - 1,3 cm, 1,9 cm and 2,5 cm widths.
500P-LAN	Can hold 15 rolls of Anesthesia tape or label products. Includes inserts, hanging bars and bottom suction cups for secure placement. Size: (H) 8,9 cm x (W) 16,3 cm x (D) 11,4 cm <i>Meant for tape and label products</i>
500P-LAN-LBL	Can hold 15 rolls of Anesthesia label products. Includes inserts, hanging bars and bottom suction cups for secure placement. Size: (H) 8,9 cm x (W) 16,3 cm x (D) 11,4 cm <i>Meant for label products only</i>