

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Spinfin® Double Stir Bars
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2024-01-03

VWR Catalogue Number	Length in mm [inches]	Diameter in mm [inches]
58948-998	12 mm (½")	12 mm (½")
58949-004	15 mm (⅝")	19 mm (¾")
58949-000	19 mm (¾")	25 mm (1")

Quality System Compliance

Products are manufactured under the [ISO 9001:2015](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Autoclavable at 121° C for 20 minutes.

Tolerance ± 0.020"

Product Specifications

Material	Dupont White Teflon #807N-X; Magnet, ALNICO 5; Meet FDA Regulations 21CFR177.1550, 21CFR177.1520 and Teflon used to totally encapsulate the magnetic bar meets US Pharmacopeias Class VI. This product does not contain or come in contact Formaldehyde or other contaminants.
Sterilization	The product is non-sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.

BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	The product is Latex free.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP), DEHP (di (2- ethylhexyl) phthalate)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2023-06-14 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Store at ambient warehouse temperature, low humidity, protected from direct sunlight and away from harsh chemicals.
Shelf Life	There is no definite shelf life or expiration for the product.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	January 2024