



Certificate of Conformity

VWR Bottle, Wide Neck, Square, with Screw Cap

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No: 215-1904

Lot No: VXP280819

Capacity (ml)	Thread	WxDxH (mm)	Type	Pk
4000	100-415	147x147x290	Square	1

Material: Bottle: HDPE, 100% Virgin FDA 21 CFR Compliance
Closure: PP, 100% Virgin, FDA 21 CFR Compliance

Date of Manufacture: 28th August 2019

Quality System Compliance

Products are manufactured, inspected and found to be in compliance with product and quality specification requirements as documented in the ISO 9001 :2015 & ISO 13485 :2016 Quality Management System. Inspection records are reviewed and signed for product release.

QC Testing

Representative production samples are collected and inspected in accordance with current applicable standard protocols



Certificate of Conformity

Test Name	Test Reference	Accuracy	Result
Leak Test	QWI-QC-04	General Insp. Level-1, AQL -1.5	Leak Proof
Autoclaving	QWI-QC-07	General Insp. Level-1, AQL -1.5	NA
Dimensional	QSF-QCIR-02	General Insp. Level-1, AQL -1.5	Within Tolerance

BSE/TSE Statement: The above mentioned product is manufactured without the use of Raw materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge.

Non-Pyrogenic: Not applicable

DNase and RNase Free: Not applicable

Cytotoxicity: Not applicable

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

This document has been produced electronically and is valid without a signature

Date of issue: **October 2020**