


Recombinant Human Adenosine Receptor A2a Protein (Fc Tag) (A318302)

Specifications:

Name:	Recombinant Human Adenosine Receptor A2a Protein (Fc Tag)
Description:	This protein is constructed from the following Adenosine Receptor A2a sequences: (Met1-Ser7) (Ser67-His tag77) (Asn144-Pro173) (Cys259-Pro266). The human Fc tag consists of amino acids Glu99-Ala330.
Applications:	SDS-PAGE
Expression System:	HEK293 cells
Nature:	Recombinant
Protein Species:	Human
Protein Length:	Protein fragment.
Sequence:	ADORA2A(Met1-Ser7)(Ser67-HisTag77)(Asn144-Pro173)(Cys259-Pro266)+hFc(Glu99-Ala330)
Tag:	C-terminal Human Fc Tag
Molecular Weight:	The protein has a predicted molecular mass of 31.7 kDa after removal of the signal peptide. The apparent molecular mass of ADORA2A-hFc is approximately 35-55 kDa due to glycosylation.
Conjugate:	Unconjugated
Purity:	> 95% by SDS-PAGE and Coomassie blue staining.
Product Form:	Lyophilized
Concentration:	Reconstitution dependent.
Formulation:	Lyophilized from sterile Phosphate Buffered Saline, pH 7.4. Normally 5%-8% Trehalose is added as a protectant before lyophilization.
Storage:	Shipped at 4°C. Lyophilized: Store at -20°C to -80°C. Reconstituted: Aliquot and store at -80°C. Product is stable for one year. Avoid freeze/thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

Recombinant Human Adenosine Receptor A2a Protein (Fc Tag) (A318302)

Images:



Recombinant Human Adenosine Receptor A2a Protein (Fc Tag) (A318302) on SDS-PAGE under reducing conditions.