

	Document Number	Revision	Page	Date
	TCB002	A0	1 of 1	07-20-2014
PRODUCT DATA SHEET		2000ml Tissue Culture Bottles		

1: Material:

Bottle: Polystyrene PS Color: Clear


Cap: High Density Polyethylene Color: Purple 259

Membrane: 0.22µm Hydrophobic (Only use for Vented Cap)

2: Related Catalog Numbers:


VWR EU EAN. No.	VWR NA Cat. No.	Capacity	Cell growth area	Cap Style	Sterile	Surface Type	Packout Description
734-2753	10769-802	2000ml	1900 cm ²	Plug seal	Yes	Treated	5/ Bag, 40/Case

3: Dimensions=mm




Printing color: Purple 259

plug seal cap



vented cap



4: Features

- Temperature range: -20°C to +50°C
- Autoclavability: No
- Non-pyrogenic
- Every package case has labeled lot No.
- Every bottle has silk-printed graduation
- Shelf life: 3 years after month of production

5: STERILITY

- Sterilized by gamma irradiation SAL10⁻⁶ (ISO11137)

6: PRODUCT ISO 9000 CERTIFICATION

- Manufactured in a class 100,000 room environment.
- Manufactured under EN ISO13485 :2003 and ISO9001:2008 quality management system.

Drawn	Approved
Date 10 July 2014	Date 10 July 2014
Name Alex Fang	Name James