

Certificate of Analysis

Product: STERILE WATER FOR IRRIGATION, USP - 500 ML

Lot #: G109504
Code: 2F7113
Expiry: OCT 2016

Manufactured Date: 10/16/2013
Sterilization Date: 10/16/2013

TEST NUMBER	TEST NAME	LIMIT	RESULT
11-20-15-001	Oxidizable Substances	Pass	Pass
11-21-16-003	pH at 25 deg. C	5.0 - 7.0	5.8
11-22-03-007	Conductivity Determination of Water (USP/EP)	NMT 5 uS/cm	1 uS/cm
20-03-07-006	Particle Analysis - Large Volume (Instrumental)		
	Particles $\geq$ 10 micron	NMT 25 particles/mL $\geq$ 10 microns	1 ct/mL
	Particles $\geq$ 25 micron	NMT 3 particles/mL $\geq$ 25 microns	0 ct/mL
13-01-V	Bacterial Endotoxin	< 0.25 EU/mL = Pass	Pass
14-04-01-120	Sterility	Pass Parametric Release	Pass
	PARAMETRIC RELEASE IN ACCORDANCE WITH U.S. FDA REQUIREMENTS		

THIS PRODUCT WAS MANUFACTURED ACCORDING TO CURRENT GMP,
FULLFILLS THE REQUIREMENTS OF THE MASTER PRODUCTION DOCUMENT AND TESTED TO USP COMPLIANCE.

E. Soliz 10-23-13
Prepared By/Date

ELLEN SOLIZ

TPoutter 10-24-13
Quality Management Approval /Date

TPoutter