

Material	23367.321
Material description	Dichloromethane
Grade	GPR RECTAPUR Stabilized with about 0,002 % of 2-methyl-2-butene
CAS Number	75-09-2
Molecular formula	CH ₂ Cl ₂
Molecular mass	84.93

Characteristics	Specifications
Assay	≥ 99 %
Stabilizer (Isoamylene)	Confirmed
IR Spectrum	Passes test
Boiling point	39.0 - 40.5 °C
Acidity	≤ 0.0004 meq/g
Density (20/4)	1.322 - 1.327
Evaporation residue	≤ 10 ppm
Free chlorine	≤ 0.3 ppm
Heavy metals (as Pb)	≤ 2 ppm
Water	≤ 0.02 %
Cl (Chloride)	≤ 3 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.