

COWIE TECHNOLOGY GROUP LTD
RIDGEWAY, COULBY NEWHAM
MIDDLESBROUGH, CLEVELAND
TS8 0TQ, UK

Tel: +44 (0)1642 599190

technical@cowie.com
www.cowie.com

Item Description:	Octahedral Stirrer Bar 76mmL x 13mmDiam
Item Reference Number:	FT.8360
Customer:	VWR International BVBA
Order No:	4513949075
Date:	23 rd November 2021
Attention: Marco Wohlschlegel	Email: Marco.wohlschlegel@avantorsciences.com

We confirm that the item listed above complies with the Item Code/Reference specification.

Mfg Job No: J110025 **Quantity:** 200 **Date of Manufacture:** November 2021

Material: TF1620 **Manufacturer:** 3M™ Dyneon™ **Batch No:** 7973788

Ref. No.	Characteristic	Spec. / Tolerance	Test Method	Quantity / Frequency
1	Surface Finish	No Visual Defects	Visual Standard	25 Parts / Every Batch
2	Appearance	No Dirt No Chips No Cracks No Flash No 'Chalky' Areas	Visual Standard	All Parts / Every Batch
3	Length	72 ~ 80 mm	Visual/Vernier	25 Parts / Every Batch
4	Width	11.7 ~ 14.3 mm	Visual/Vernier	25 Parts / Every Batch
5	Height	11.7 ~ 14.3 mm	Visual/Vernier	25 Parts / Every Batch
6	Batch Quantity	Correct to Order	Visual	Every Batch
7	Porosity	No Visible Dye	Dye Penetration	1 Part / Batch
8	Magnetism	Support Weight of Stirrer Bar	Suspension against Ferrous Plate	1 Parts / Bag

THIS DOCUMENT IS VALID ONLY TO THE CUSTOMER IDENTIFIED ABOVE

Cowie Technology Group Ltd designs & manufactures under ISO9001:2015 certified by Lloyds Register QA.
PTFE used in the manufacture of this product is VIRGIN PTFE GRADE TF1620 manufactured by 3M™ DYNEON™

3M™ DYNEON™ TF1620 – COMPLIANCE STATUS

FDA – Directive FDA 21 CFR 177.1550

USP Class VI – USP 27, NF 22, 2004 for Plastics at 70°C.

RoHS 2 – Directive 2011/65/EU & Dangerous Substances – 76/769/EC

RoHS 2 Phthalates – Directive 2011/65/EU as amended by 2015/863/EU

EU REACH Regulations – SVHC Current Candidate List

Bisphenol A – Absent

Latex – Absent

ADCF – Animal Derived Components Free, BSE/TSE – EMEA/410/01/rev.3, GMO Free

Endotoxins – Free

POP, PFOS, PFOA including Homologues, Derivatives and Salts of PFOA, APFO, DEHP, Phthalates, VOC Free

ODC – Class 1


Dodd-Frank Act – Gold, Coltan, Cassiterite, Wolframite and their derivatives originating from the Democratic

Republic of Congo – Absent

California Proposition 65 – Absent of chemicals contained on the California Proposition 65 list

It is the responsibility of the user to ensure that the products are suitable for their intended use

Signed


QUALITY MANAGEMENT REPRESENTATIVE