



VWR® N95 DISPOSABLE RESPIRATORS

As per the manufacturer, the below product meets the following criteria:

VWR Catalog Number:	89201-508
Lot Number:	022218
Description:	VWR N95 RESPIRATOR PK20
Date of Manufacture:	Feb 22, 2018
Date of Expiration:	Feb 21, 2021
Country of Origin:	China

Quality System Compliance

Products are manufactured under the ISO 9001:2008 & ISO13485:2003 standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP. Inspection records are reviewed and signed off by qualified personnel for product release.

QC Testing

Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E and analysis in accordance with IEST-RP-CC004.3

PFE **Spec:** 95% **Test Result:** 99.714% ~ 99.869%(Pass)

Delta P **Spec:** 35mmH₂O **Test Result:** 7.8 ~ 8.9 mmH₂O(Pass)

Product Specifications

Material:	100gPE/38gMB/160gPE
BSE/TSE Statement:	The above product contains no ingredients of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Spongiform Encephalopathy (TSE)/Bovine spongiform Encephalopathy (BSE).
Non-Pyrogenic:	N/A
DNase & RNase Free:	N/A
Cytotoxicity:	N/A
Latex Free Statement:	The above product is latex free.

Signed:



Jamie Ethier
VP Global Quality
VWR, Part of Avantor

Date: October 12, 2020